




APIRT02

# Infrared Forehead Thermometer

## Instruction manual



Please read the manual carefully before you use the unit, and keep for future reference.

### Introduction

▲Your new Infrared forehead thermometer uses advanced infrared (IR) technology to measure temperature instantly and accurately on the forehead or object.

▲Easy to use and less measurement time

This thermometer safely and hygienically measures a persons temperature without contact to the body. Its ergonomic design makes this thermometer, simple and easy to check a persons temperature. It only takes 1 second to take an accurate measurement or reading.

▲Body and Object modes

This thermometer is able to measure temperature of a body or object. Measuring range of object mode is 0.0 to 100.0°C (32.0°F~199.9°F ). That means apart from to taking a persons body temperature, it also can measure the surface temperature of object. Examples below

Surface temperature of milk in a babies bottle;

Surface temperature of a babies bath;

Environment temperature.

▲LED colour indicator and alarm sounds

When body temperature is over 37.5°C, the colour indicator will show red and an alarm sounds.

▲Memory function

It stores up to 20 sets of recent measurement data.

▲Auto power off

Automatically power off if left idle for Approx 60s.

▲Measurement

Measurement time interval is 1 second and the measurement distance is within 1~3cm.

▲Intended use

The Infrared Thermometer is intended for the intermittent measurement and monitoring of human body temperature from the forehead or the surface temperature of an object. The device is intended for use by people of all ages, for homecare or in hospitals.

### Safety Information

To ensure the correct use of this product, all basic safety measures should always be followed including the warnings and cautions listed in the instruction manual.

Symbol descriptions
The following symbols may appear in this manual, on the label, on the device, or on it's accessories. Some of the symbols represent standards and compliances associated with the device and its use.

	Regulatory compliance logo for Australia and New Zealand. Covers Electrical Safety and EMC.
	WARNING: This alert identifies hazards that may cause serious personal injury or death.
	CAUTION: This alert identifies hazards that may cause minor personal injury, product damage, or property damage.
	Type BF applied part
	Manufacturer
<b>SN</b>	Specifies serial number
	<b>DISPOSAL:</b> Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.
	Direct current
	Follow instructions for use
	CAUTION: Consult accompanying documents

▲This thermometer is not intended to substitute for a consultation with your physician. The forehead scan temperature serves as a reference only. It cannot be a judgment on fever.

▲Basic safety precautions should always be observed, especially when the thermometer is used on or near children and disabled persons.

▲Please place the device out of reach of children.

▲Avoid direct sunlight.

▲Do not touch the lens.

▲No modification of this device is allowed.

▲The swallowing of small parts like packing bag, battery, battery cover and so on may cause suffocation.

▲Please do not use a dilution agent, alcohol or petrol to clean the unit.

▲Please treat device gently and prevent it from falling from a high place.

▲Please do not immerse it in liquid.

▲Never leave the batteries in the battery compartment for an extended time without use, as they may leak and cause damage to the unit.

▲Please take out the batteries if you do not intend to use within 3 months.

▲Replace with new batteries if the unit displays a low battery symbol.

▲Do not mix the old and new batteries together.

▲ **WARNING:**

Do not dispose of electrical appliances as unsorted municipal waste, use separate collection facilities. Contact your local government for information regarding the collection systems available. If electrical appliances are disposed of in landfills or dumps, hazardous substances can leak into the groundwater and get into the food chain, damaging your health and well-being.


### Classification

1. Internally powered equipment;
2. Type BF applied part; .
3. Protection against ingress of water or partial matter: IP21;
4. Not category AP/APG equipment;
5. Mode of operation: Continuous operation.

▲The user must check that the equipment functions safely and that it is in proper working condition before being used.


### Product Structure

#### -Display



1. Body temperature indicator
2. Object temperature indicator
3. Memory mode
4. Temperature reading
5. Temperature unit
6. Low battery indicator

#### -Body



1. Auto-sense Probe
2. Display LED
3. Display Screen
4. SET Button
5. Memory Button
6. START/STOP Button (Power ON/OFF Measure)
7. Battery Cover

### Battery Installation

Remove the battery cover from the battery compartment, insert the battery.

1. Remove the battery cover from the battery compartment as the arrow direction accordingly.



2. Insert 2 AAA Alkaline batteries into the compartment and ensure each battery is in the proper direction, Positive (+) and Negative (-) are displayed on the back of battery cover.



3. Close the battery cover.



### Low battery and replacement

If the low battery symbol appears on the screen you must replace with new batteries before the unit will work again.

### Battery type and replacement

Please use 2pcs AAA identical 1.5V Alkaline batteries.

Do not use the batteries beyond their expiry date.

Please remove the batteries if you do not need the unit for an extended period.

### WARNING


Dispose of the battery in accordance with all federal, state and local laws. To avoid fire and explosion hazard, do not burn or incinerate the battery.

### Setting Mode


#### How To Set

##### 1. Mode setting:

When powered press SET button (after the unit beeps) and the screen will display or . Mode will change between and when the SET button is pressed. It will display the Body Mode automatically when you switch to , and display the Object Mode automatically when you switch to .



Body mode




Object mode


Note: The Body Mode is used to measure the forehead temperature. Object Mode is used to measure the object temperature.

##### 2. Temperature Unit setting:

With the power on, long press the SET button for 5 seconds (approx) until you hear 1 beep. The screen will display or only. Press the SET button again to change between Celsius & Fahrenheit. Then press power button twice to confirm. The unit is then ready to measure temperature



Celsius/°C



Fahrenheit/°F

### Using your device

#### Pre-measurement

About Normal Body Temperature & Fever

The temperature in the forehead and temple area differs from the internal temperature, which is taken orally or rectally. Vasoconstriction, an effect which constricts the blood vessels and cools the skin, can occur during the early stages of a fever. In this case, the temperature measured by the Infrared thermometer may be unusually low. If the measurement therefore does not match the patient's own perception or is unusually low, repeat the measurement every 15 minutes. As a reference, you can also measure the internal body temperature using a conventional oral or rectal thermometer.

Body temperature can vary from one individual/person to next.

It also varies by location on the body and time of day. Below shows the statistical normal ranges from different sites.

Please keep in mind that temperatures measured from different sites, even at the same time, should not be directly compared. Fever indicates that the body temperature is higher than normal. This symptom may be caused by infection, overdressing or immunization. Some people may not experience fever even when they are ill.

These include, but are not limited to, infants younger than 3 months old, persons with compromised immune systems, persons taking antibiotics, steroids or antipyretics (aspirin, ibuprofen, acetaminophen), or persons with certain chronic illnesses. Please consult your physician when you feel ill even if you do not have fever.

#### Table\*1 Body Site Normal Temperature Range

Body Site	Normal Temperature Range
Oral	0.6°C (1°F) or more above or below 37°C (98.6°F)
Rectal/ear	0.3°C to 0.6°C (0.5°F to 1°F) higher than oral temperature
Axillary (armpit)	0.3°C to 0.6°C (0.5°F to 1°F) lower than oral temperature

**Note:** Body Temperature at WebMD;

**website:** <http://firstaid.webmd.com/body-temperature>; retrieved at 2010 Jan 7.


#### As a Body Thermometer

1). Press the START/STOP button, all symbols appear on the display, you can hear 1 short beep, then unit °C/°F on the display start to flash.

2). Press the SET button to select Body Mode, the temperature unit flashes.

3). Move the probe close to the forehead, make sure the probe is flat and take measurement with a distance between 1~3cm. Press the START/STOP button, the meter will take measurement.

4). The measurement will be finished in 1 second, when it has been completed you can hear 1 short beep, in the meantime, the reading will appear on the display with green LED blinks.



**Note:**

● If the reading is < 37.5°C (99.5°F), the display will show together with green LED.

● If the reading is ≥ 37.5°C (99.5°F) and < 43°C (109.4°F), the display will show together with red LED and ten short beeps.

● Only after the unit °C/°F start to flash again, it can continue to take measurements.

● As the forehead measurement temperature is likely to be affected by sweat, oil and the surroundings, the reading shall be taken as a reference only.


● If the probe is placed at an angle close to the forehead measurement, the reading will be affected by surrounding temperature.

● Babies' skin reacts very quickly in the ambient temperature. Therefore, do not take their temperature with the non-contact thermometer during/after breastfeeding, because the skin temperature may be lower than the internal body temperature.



**As a Object Thermometer**

- 1). Press the START/STOP button, all symbols appear on the display , you can hear 1 short beep, then unit °C/°F on the display starts to flash.
- 2). Press the SET button to select Object Mode, the temperature unit flashes.
- 3). Move the probe close to the object, make sure the probe is flat and take measurement with a distance between 1~3cm. Press the START/STOP button, The thermometer will take a measurement.
- 4). The measurement will be finished in 1 second, when it has been completed, you can hear 1 short beep, in the meantime, the reading will appear on the display with a green LED.



**Note:**  
 ●As the object scan temperature is likely to be affected by the surroundings, the reading shall be taken as a reference only.

- If the probe is placed at an angle close the object scan, the reading will be affected by surrounding temperature.
- Take care "Polished metals like a stainless steel kettle or saucepan will reflect the surrounding temperature of the room rather than their own emitted temperature. You may get a reading of 38°C (100°F) on the side of a pot on the stove and think it is safe to touch but you will burn your fingers as it could be boiling at 100° (212°F) The same is true for very cold materials as well so please take care when using this device on reflective surfaces"

**Memory-recall of measurements**


This Non-contact forehead thermometer automatically stores 20 sets measurement value, the oldest record will be replaced by the latest measurement value when more than 20 sets is exceeded.

**Recall memory**

Press the MEM Button to enter the memory mode when power off. Each time you press the MEM button, a number (from 1 to 20) will be displayed along with the symbol M, 1 second later, the measurement will be shown.

Subsequent measurements can be displayed, one after the other by pressing the MEM Button each time.

**Note:**  
 The reading number "1" is the latest measurement, "20" is the oldest measurement. When the oldest measurement is read, pressing the MEM Button again will read the reading "1" again.



**Memory-Clear Of Measurements**

If you want to remove all stored memories permanently, Long press the MEM Button for 10 seconds, until CLR appears. CLR will flash and the unit will beep 5 times and then switch off.



**Exceptional Situation**

Symbol	Explanation
	In Body Mode, Measured temperature is above Measuring range of 43°C/109.4°F it is shown on the LCD with red LED.
	In Body Mode, Measured temperature is below the Measuring range of 32°C/89.6°F, it is shown on the LCD with green LED
	In Object Mode, Measured temperature is above Measuring range 100.0°C/199.9°F or environmental temperature is above system Operating range 40°C/104.0°F, it is shown on the LCD with red LED.
	In Object Mode, Measured temperature is below Measuring range 0.0°C/32.0°F, or environmental temperature is be low system operating range 5°C/41.0°F, it is shown on the LCD with green LED.

	Low battery, replace all the flat batteries with new ones.
	Thermometer system fails or affected by electric magnetic field

⚠ Please contact the distributor if you can not solve the problem, do not disassemble the unit by yourself!

**Care And Maintenance**

- Care for the main unit**
- Keep the unit in the packaging box when not in use.
  - Clean the unit with a soft dry cloth. Do not use any abrasive or volatile cleaners.
  - Never immerse the unit or any of it's components in water.

**Maintenance**

Do not clean the unit with chemicals, thinners or gasoline etc.	Store the unit in a clean and dry location. Do not subject the unit to extreme hot or cold temperatures, humidity or direct sunlight.
Remove the batteries if the unit will not be used for 3 months or longer.	Do not use the unit under the influence of electromagnetic interference (nearby cellphones, microwave etc.)

**Note:** We will not be responsible for any quality problem if you do not care and maintain the product as instructed.

**Specifications**

Description	Non-contact Infrared Forehead Thermometer
Display	LCD digital display
Measuring localization	Forehead and object surface
Measurement range	Body mode 32.0°C~43.0°C (89.6°F~109.4°F); Object mode 0.0~100.0°C (32.0°F~199.9°F);
Temperature unit	°C/°F
Display resolution	0.1°C/°F
Accuracy	±0.2°C/±0.4°F (within 36.0°C~39.0°C/ 96.8°F~102.2°F)
Memory function	20 sets memory of measurement values
Beep alarm	One short beep when power on and start measurement One long beep when measurement reading is below 37.5°C/99.5°F. 10 short beeps when measurement reading is greater than or equal to 37.5°C/99.5°F 3 short beeps when system fails
LED color indicator	Green: Temperature < 37.5°C/99.5°F Red: Temperature ≥ 37.5°C
Power source	2 x AAA Alkaline batteries 1.5V (batteries not included)
Automatic power-off	In 60s Manual shutdown In 8s
Main unit weight	Approx. 75g (without batteries)
Main unit size	L160 mm x W55 mm x H38 mm
Battery life	Could be used for 300 times for normal condition (Alkaline batteries)
Accessories	Instruction manual
Operating environment	Body mode 10.0°C~40.0 °C/50.0°F~104.0 °F Object mode 5°C~40.0°C/41.0°F~104.0°F
Storage environment	Temperature: -20°C~+50°C/-4.0°F~+122.0 °F Humidity: 15% ~95% avoid crash, sun burn or rain during transportation

**Clinical Measurement Accuracy And Safety Verification:**

The product has passed clinical trials. The measured results of the infrared forehead thermometer was compared with the measured results of mercury thermometers, the deviation average  $\Delta t_b = 0.011^\circ\text{C}$  not exceeding 0.3°C; the clinical repeatability of the infrared forehead thermometer SR = 0.100°C, not exceeding 0.3°C. The measured results up to the laboratory standard and the clinical standard. Therefore, the deviation average and the clinical repeatability of the infrared forehead thermometer complied with the regulatory requirement ISO 80601-2-56. The conclusions are drawn from the clinical trials, the accuracy and safety complied with the regulatory requirement.

**Safety information & warnings**

- ▲ Do not use the unit in damp working conditions; the unit is not waterproof.
- ▲ Do not try opening the shell or trying to repair the unit by yourself. If the unit has any problems, please contact or take it to an authorised repair agent for repair.
- ▲ Do not use a chemical substance, cleaning solvent or cleaner to clean the unit, please use a moist soft cloth.

**EMC**

**Table 1**  
Guidance and manufacturer's declaration-electromagnetic emissions

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Infrared forehead thermometer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class [B]	The Infrared forehead thermometer is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

**Table 2**  
Guidance and manufacturer's declaration-electromagnetic emissions

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV input/output lines: ±1 kV	Power supply lines: ±2 kV input/output lines: ±1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	line(s) to line(s): ±1 kV. line(s) to earth: ±2 kV. 100 kHz repetition frequency	line(s) to line(s): ±1 kV. line(s) to earth: ±2 kV. 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 300 cycle	0% 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 300 cycle	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U<sub>r</sub> is the a.c. mains voltage prior to application of the test level.

**Table 3**  
Guidance and manufacturer's declaration-electromagnetic emissions

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC61000-4-6	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the Infrared forehead thermometer, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distances: $d=0.35\sqrt{P}$ , $d=1.2\sqrt{P}$ , $d=2.3\sqrt{P}$
Radiated RF IEC61000-4-3	10V/m, 80% Am at 1kHz	10V/m, 80% Am at 1kHz	80MHz to 800MHz: $d=1.2\sqrt{P}$ , 800MHz to 2.7GHz: $d=2.3\sqrt{P}$ ; Where, P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, d is the recommended separation distance in meters (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.  
 NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

▲ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Infrared forehead thermometer is used exceeds the applicable RF compliance level above, the Infrared forehead thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Infrared forehead thermometer.  
 ▲ Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Table 4**  
Recommended separation distances between portable and mobile RF communications equipment and the

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=3.5\sqrt{P}$	80MHz to 800MHz $d=1.2\sqrt{P}$	800MHz to 2.7GHz $d=2.3\sqrt{P}$
0.01	/	0.12	0.23
0.1	/	0.38	0.73
1	/	1.2	2.3
10	/	3.8	7.3
100	/	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.  
 NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.  
 NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**Table 5**  
Guidance and manufacturer's declaration-electromagnetic emissions

Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation n b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380-390	TETRA 400	Pulse modulation b) 18 Hz	1.8	0.3	27
450	380-390	GMSR 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation b) 217 Hz	0.2	0.3	9
745						
780						
810	800-960	TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25, UMTS	Pulse modulation b) 217 Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0.2	0.3	9
5240						
5785						

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.  
 b) The carrier wave is modulated using a 50% duty cycle square wave signal.  
 c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVEL .S shall be calculated using the following equation:  
 $E = \frac{E_{ref}}{d}$   
 Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

**Support Information/Authorized agent**

Imported by TDJ Australia PTY LTD. ARTG Ref No: 333301

Please retain your proof of purchase receipt in order to make a claim about faulty products or services. If you need assistance setting up or using your Aerpro product, or make a claim now or in the future, call Aerpro Support. Australia

TEL: 03 – 8587 8898

FAX: 03 – 8587 8866

Mon-Fri 9am – 5pm AEST

If you would like to download a digital copy of this manual, or other Aerpro manuals/software, please visit the <http://aerpro.com> website and search for the product model number for more information, accessories and products.

This manual is considered correct at time of printing but is subject to change. For latest manuals and updates refer to the website.

Model:UFR106

**Manufacturer**  
 Shenzhen Urion Technology Co.,Ltd.  
 Floor 4-6th of Building D,Jiale Science&Technology Industrial Zone,No.3,ChuangWei Road,Heshuikou Community,MaTian Street, GuangMing New District,518106 ShenZhen,PEOPLE'S REPUBLIC OF CHINA  
 Tel: (86) -755-29231308 E-mail:urion@urion.com.cn  
 MADE IN CHINA

Shanghai International Holding Corp. GmbH (Europe)  
 Eiffestrasse 80, 20537 Hamburg, Germany  
 Tel:+49-40-2513175 E-mail:shholding@hotmail.com

